



Making Waves in Medical Therapy!

Alpha Waves

Volume I, Issue 2

December 2005

EXPANDING THE KNOWLEDGE ABOUT ALPHA-STIM® FOR ANIMALS

The Wild West Veterinary Conference held Oct 5-9, 2005 in Reno, Nevada was yet another milestone for Alpha-Stim in the veterinary market with record sales for EPII! And that same month *Equine Veterinary Science Journal* published "Microcurrent Electrical Therapy" written by Dr. Ava Frick with technical application by Doug McCauley on the use of MET treatment of a recalcitrant wound in a mare. The article is available on the website.

"We are building a great educational format for doctors, distributors, and clients using Alpha-Stim® for animal applications. With research, seminars and an upcoming handbook on animal protocols, everyone can be utilizing this fantastic therapeutic methodology with certainty." Says Dr. Ava Frick who according to Dr. Dan Kirsch of EPII is the world's leading authority of microcurrent therapy in animals.

March 31 and April 1 in Union, Missouri Dr. Ava Frick will instruct veteri-

narians, doctors, nurses, technicians, and Alpha-Stim® distributors on the theories and application of Cranial Electrotherapy Stimulation (CES) and Microcurrent Electrotherapy (MET). The program is 11 hours; 8 of lecture and 3 of hands-on usage on your choice of dogs, horses, and humans. This seminar is approved for 11 CE units.

This is a great opportunity to learn about real cases, real results, and how to get the most out of your time and equipment. Added entertainment is scheduled for Friday evening and Saturday afternoon or plan a little extra time in Missouri to visit some of the local wineries.

Early bird registration closes March 15, 2006 so call now to get signed up!

Call either 636-583-1700 or 866-836-3900 to get fees, lodging, and to register. See you there!

WHAT SETTINGS DO YOU USE ON DOGS?

Animals respond to microcurrent therapy much better than most humans. Younger bodies also need less intensity and duration than older ones. My standard protocol for adult dogs with CES is 0.5 Hz and 100 to 300 microamperes. For anxiety and OCD disorders start at 1 hour per treatment reducing to 20 minutes as they respond. In pain conditions follow MET with 20 minutes of CES (dogs average 3 minutes post relaxation at 9-15 minutes). MET treatment with probes for pain is set at 400 to 600 microamperes as the patient will tolerate.

Upcoming Alpha Waves ... Dr. Frick having directed more than **3,000 animal Alpha-Stim® treatments will continue to address specific protocols to customize each species you work with!**

What people are saying about Dr. Frick's seminars.

"Dr. Frick did an excellent job at presenting, explaining, giving examples and visuals both hands-on and in notes. I think it would be good to make more doctors aware of Alpha-Stim therapy as an alternative to drugs."

Sheila Snodgrass, DVM

"Dr. Frick is not only expert, but a very good teacher, too."

Jerry Fila, Jr., MBA

"She was great - the best vet seminar I have attended"

Dyanna McCain, DVM

REMEMBER:

When treating animals rewet the probe pads every minute with Alpha Conducting Solution.

Making Waves in Medical Therapy!

1825 Denmark Road
Union, MO 63084

Phone: 636-583-0094
1-888-583-0944
Fax: 636-584-0371
E-mail:

dravafrick@charterinternet.com



STOP BY OUR EXHIBIT AT THE MISSOURI VETERINARY MEDICAL ASSOCIATION CONFERENCE AT BRANSON JANUARY 19-20, 2006 AND GET A FREE TREATMENT ON YOURSELF FOR WHATEVER PAIN AILS YOU!

When the fur goes, where do you go for help?

History: Kate is a four year old cat who up until 18 months ago was the perfect house pet. The problems began when the owners installed bamboo flooring into their home. For days Katie had to listen to the repetitive sound of the electric hammers nailing in the floor. Shortly thereafter the owners noticed a balding pattern to her abdomen, tail and rear legs. Eight months later they realized that she was also urinating behind the piano. Needless to say, this brought them and her to the veterinarian. Four visits and injections of corticosteroids failed to resolve the etiology of her problem, which was anxiety and obsessive compulsive disorder.

Katie before Alpha-Stim® therapy:



Katie after 60 days of Alpha-Stim® therapy:



Katie's Alpha-Stim® therapy protocol:

1 hour CES per day for 2 weeks then 20 minutes per day for 2 weeks then every other day. On reduction of the frequency she began licking herself more so the treatments were increased to daily administration. Gradually over time the interval was again extended. Abnormal urination habits ceased and fur returned.